


# TID PHAROS

## ILLUMINATED DEPRESSOR




- TID Pharos is the market's ONLY one-piece, complete illuminated depressor that connects directly into your vitrectomy machine
- Conventional, ball-point depressor tip with a "Side Fire Illumination" feature allows for a better view, less glare and improved surgeon ergonomics
- The TID Pharos is designed to offer the surgeon the ability to simultaneously perform unassisted scleral depression while shaving vitreous base
- Additional Benefits include: identifying retained lens fragments, identifying radiation seed placements and locating breaks/tears in the outer periphery



# TID | TODORICH

## ILLUMINATED DEPRESSOR

- First generation, illuminated depressor cap-over design that slides directly onto your vitrectomy machine endoilluminator
- Perform unassisted scleral depression while shaving the vitreous base
- Compatible with Alcon, Bausch & Lomb and DORC EVA standard endoilluminators



Catalog Numbering and Ordering information are located on the reverse side





# Catalog Numbers and Descriptions

## **TID PHAROS ILLUMINATED DEPRESSOR**

- VS0801A TID Pharos Illuminated Depressor - Alcon, Sterile, 5/bx
- VS0801B TID Pharos Illuminated Depressor - Bausch and Lomb, Sterile, 5/bx
- VS0801D TID Pharos Illuminated Depressor - DORC EVA, Sterile, 5/bx

**Developed in cooperation with Bozho Todorich, MD**

## **TID ILLUMINATED DEPRESSOR CAP**

- VS0800 TID Illuminated Depressor Cap, Sterile, 5/bx

**Developed in cooperation with Bozho Todorich, MD**

## **ORDERING INFORMATION:**



**TID/TID PHAROS**

Vortex Surgical  
4 Research Park Drive, Suite 124  
Saint Charles, MO 63304 USA

Phone: 636-778-4350  
Email: [info@vortexsurgical.com](mailto:info@vortexsurgical.com)  
Web: [www.vortexsurgical.com](http://www.vortexsurgical.com)



**VORTEX SURGICAL**