

INJECTION & DRAINAGE

Vortex Surgical offers a wide range of unique injection and drainage devices to fit any surgeon's preference and every patient's needs

EXTERNAL DRAINAGE AND DEPRESSION - EDD



- 28ga drainage needle x 2.4mm long
- Eliminates the need for a cut down while buckling
- Precise localization of drainage point
- Fixed drainage needle length for improved control
- Helps reduce the risk of retinal incarceration and hemorrhages
- Straight and Curved versions available

NANO SUBRET INJECTION DEVICE

- Optimized for use without vitrectomy
- Transcleral / transvitreal passage
- Precise subretinal delivery of medications, gene therapy and stem cells
- Ideal for tumor biopsies
- Useful in removing retained subretinal PFO
- 28ga self-penetrating outer needle with 41ga beveled injection needle
- Straight and Curved versions available



INTERNAL DELIVERY DEVICE - IDD

- 25ga shaft with 27ga extendable nitinol
- Beveled end for easy penetration of pathology
- Ideal for delivery of viscoelastics
- Micro-does tubing set included to minimize injection waste.
- Available in 25ga



Catalog Numbering and Ordering information are located on the reverse side





EXTERNAL DRAINAGE AND DEPRESSION

- VS0290 External Drainage and Depression Device - Straight, Sterile, 5/bx
VS0290C External Drainage and Depression Device - Curved, Sterile, 5/bx

Developed in cooperation with **Wilken Parke MD, Emmanuel Chang MD, Gholam Peyman MD and Charalampos Liver-Rallatos MD**

NANO SUBRET INJECTION DEVICE

- VS0225 NANO SubRet Injection Device - Straight, Sterile, 5/bx
VS0225C NANO SubRet Injection Device - Curved, Sterile, 5/bx

Developed in cooperation with **Tamer Mahmoud MD, Edward Wood MD and Prethy Rao MD**

INTERNAL INJECTION DEVICE

- VS0220.25 25ga Internal Delivery Device, Sterile, 5/bx

ORDERING INFORMATION:



INJECTION & DRAINAGE

Vortex Surgical
4 Research Park Drive, Suite 124
Saint Charles, MO 63304 USA

Phone: 636-778-4350
Email: info@vortextsurgical.com
Web: www.vortextsurgical.com



VORTEX SURGICAL